

SAFETY DATA SHEET

1. Identification

1. Identification			
Product identifier	Prednisone Tablets		
Other means of identification			
Catalog number	1559505		
Recommended use	Specified quality tests and a		
Recommended restrictions	Not for use as a drug. Not for	or administration to	o humans or animals.
Manufacturer/Importer/Supplier/	Distributor information		
Manufacturer			
Company name Address	U. S. Pharmacopeia 12601 Twinbrook Parkway Rockville MD 20852-1790 United States		
Telephone	RS Technical Services	301-816-8129	
Website E-mail	www.usp.org RSTECH@usp.org		
Emergency phone number	CHEMTREC within US & Canada	1-800-424-9300	
	CHEMTREC outside US & Canada	+1 703-527-388	7
2. Hazard(s) identification			
Physical hazards	Not classified.		
Health hazards	Serious eye damage/eye irri	itation	Category 2B
	Reproductive toxicity		Category 2
	Specific target organ toxicity exposure	v, repeated	Category 1 (endocrine system)
Environmental hazards	Not classified.		
OSHA defined hazards	Not classified.		
Label elements			
Signal word	Danger		
Hazard statement	Causes eye irritation. Suspected of damaging fertility or the unborn child. Causes damage to organs (endocrine system) through prolonged or repeated exposure.		
Precautionary statement			
Prevention	Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Do not breathe dust. Wash thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection.		
Response	If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. If exposed or concerned: Get medical advice/attention. Get medical advice/attention if you feel unwell.		
Storage	Store locked up.		
Disposal	Dispose of contents/container in accordance with local/regional/national/international regulations.		
Hazard(s) not otherwise classified (HNOC)	None known.		
Supplemental information	Potent pharmacologically ac	tive material.	

3. Composition/information on ingredients

Mixture

Chemical name	Common name and synonyms	CAS number	%
Microcrystalline Cellulose		9004-34-6	<55
Dibasic Calcium Phosphat	e	7757-93-9	<45
Prednisone	1-Dehydrocortisone	53-03-2	<5
Magnesium Stearate		557-04-0	<2
Sodium Starch Glycolate		9063-38-1	<2
Stearic Acid		57-11-4	<2
Other components below reportable levels			0.32

4. First-aid measures

Inhalation	If breathing is difficult, remove to fresh air and keep at rest in a position comfortable for breathing. Call a physician if symptoms develop or persist.
Skin contact	Rinse skin with water/shower. Get medical attention if irritation develops and persists.
Eye contact	Rinse with water. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical attention if irritation develops and persists.
Ingestion	Rinse mouth. If ingestion of a large amount does occur, call a poison control center immediately.
Most important symptoms/effects, acute and delayed	Endocrine system effects. Potent pharmacologically active material Occupational exposure to small amounts may cause physiological effects.
Indication of immediate medical attention and special treatment needed	Provide general supportive measures and treat symptomatically. Treatment of corticosteroid overdose may include the following: Toxicity is low after acute ingestion. Gastrointestinal decontamination is generally not necessary.
General information	Remove from exposure. Remove contaminated clothing. For treatment advice, seek guidance from an occupational health physician or other licensed health-care provider familiar with workplace chemical exposures. In the United States, the national poison control center phone number is 1-800-222-1222. If person is not breathing, give artificial respiration. If breathing is difficult, give oxygen if available. Persons developing serious hypersensitivity (anaphylactic) reactions must receive immediate medical attention.
5. Fire-fighting measures	
Suitable extinguishing media	Water. Foam. Dry chemical or CO2. Use fire-extinguishing media appropriate for surrounding

	materials.
Unsuitable extinguishing media	None known.
Specific hazards arising from the chemical	No unusual fire or explosion hazards noted.
Special protective equipment and precautions for firefighters	Wear suitable protective equipment.
Fire fighting equipment/instructions	Use water spray to cool unopened containers. As with all fires, evacuate personnel to a safe area. Firefighters should use self-contained breathing equipment and protective clothing.
Specific methods	Use standard firefighting procedures and consider the hazards of other involved materials.
General fire hazards	No unusual fire or explosion hazards noted.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures	Keep unnecessary personnel away. Wear appropriate personal protective equipment. Avoid inhalation of dust from the spilled material. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing. Ensure adequate ventilation. For personal protection, see section 8 of the SDS.
Methods and materials for containment and cleaning up	For waste disposal, see section 13 of the SDS. Avoid the generation of dusts during clean-up. Sweep up or vacuum up spillage and collect in suitable container for disposal. Clean surface thoroughly to remove residual contamination.
Environmental precautions	Avoid discharge into drains, water courses or onto the ground.

7. Handling and storage

Precautions for safe handling

As a general rule, when handling USP Reference Standards, avoid all contact and inhalation of dust, mists, and/or vapors associated with the material. Clean equipment and work surfaces with suitable detergent or solvent after use. After removing gloves, wash hands and other exposed skin thoroughly. Select and use containment devices and personal protective equipment based on a risk assessment of material potency and exposure potential.

Conditions for safe storage, including any incompatibilities

Store in tight container as defined in the USP-NF. This material should be handled and stored per label instructions to ensure product integrity.

8. Exposure controls/personal protection

Occupational exposure limits

The following constituents are the only constituents of the product which have a PEL, TLV or other recommended exposure limit. At this time, the other constituents have no known exposure limits.

US. OSHA Table Z-1 Limits for Air Contaminants (29 CFR 1910.1000) Form Components Value Type Microcrystalline Cellulose PEL 5 mg/m3 Respirable fraction. (CAS 9004-34-6) 15 mg/m3 Total dust. **US. ACGIH Threshold Limit Values** Components Value Type Magnesium Stearate (CAS TWA 10 mg/m3 557-04-0) Microcrystalline Cellulose TWA 10 mg/m3 (CAS 9004-34-6) Stearic Acid (CAS 57-11-4) TWA 10 mg/m3 **US. NIOSH: Pocket Guide to Chemical Hazards** Components Form Value Type Microcrystalline Cellulose TWA 5 mg/m3 Respirable. (CAS 9004-34-6) 10 mg/m3 Total **Exposure limit values** Industrial Use Tuna Componente Value

Components	Туре	Value	
Prednisone (CAS 53-03-2)	TWA	50 micrograms/m3	
Biological limit values	No biological exposure limits noted	for the ingredient(s).	
Appropriate engineering controls	No open handling. For laboratory operations, use approved ventilation or containment system (biological safety cabinet, ventilated balance enclosure, glovebox). Control exposures to below the occupational exposure level (if available). Select and use containment devices and personal protective equipment based on a risk assessment of exposure potential. Cover all containers for solutions and slurries while being transferred.		below the mal
Individual protection measures,	such as personal protective equip	ment	
Eye/face protection	Near safety glasses with side shields, chemical splash goggles, or full face shield, if necessary. Base the choice of protection on the job activity and potential for contact with eyes or face. An emergency eye wash station should be available.		
Skin protection			
Hand protection		e or other impervious gloves if skin contact is possible ded in an organic solvent, wear gloves that provide pro	
Other	disposable sleeve covers and two p skin protection on the job activity, p	and degowning practices. Wear disposable lab coat, pair of gloves as appropriate for the task. Base the cho otential for skin contact and solvents and reagents in mmon areas (e.g., cafeterias) or out-of-doors.	oice of
Respiratory protection		tor (PAPR) with HEPA filters, disposable outerware an piratory protection appropriate to the task and the leve	
Thermal hazards	Wear appropriate thermal protective	e clothing, when necessary.	

General hygiene	Pharmacological effects may be seen with occupational exposure. Handling practices in this SDS
considerations	are recommendations for laboratory use of reference standards. Procedures for any other uses or
	quantities should be determined after an appropriate assessment.

9. Physical and chemical properties

Appearance	Appearance descriptions are general information and not specific to any USP lot.	
Physical state	Solid.	
Form	Tablet.	
Color	White.	
Odor	Odorless.	
Odor threshold	Not available.	
рН	Not available.	
Melting point/freezing point	Not available.	
Initial boiling point and boiling range	Not available.	
Flash point	Not available.	
Evaporation rate	Not available.	
Flammability (solid, gas)	Not available.	
Upper/lower flammability or expl	losive limits	
Flammability limit - lower (%)	Not available.	
Flammability limit - upper (%)	Not available.	
Explosive limit - lower (%)	Not available.	
Explosive limit - upper (%)	Not available.	
Vapor pressure	Not available.	
Vapor density	Not available.	
Relative density	Not available.	
Solubility(ies)		
Solubility (water)	Not available.	
Auto-ignition temperature	Not available.	
Decomposition temperature	Not available.	
Viscosity	Not available.	
10. Stability and reactivity		
Reactivity	The product is stable and non-reactive under normal conditions of use, storage and transport.	
Chemical stability	Material is stable under normal conditions.	
Possibility of hazardous reactions	No dangerous reaction known under conditions of normal use.	
Conditions to avoid	Contact with incompatible materials.	
Incompatible materials	Strong acids. Strong oxidizing agents. Amines. Isocyanates. Fluorine. Strong bases.	
Hazardous decomposition products	MgOx, POx, NaOx, CaOx. Irritating and/or toxic fumes or gases. Emits toxic fumes under fire conditions.	
11. Toxicological informat	ion	
Information on likely routes of e	xposure	
Inhalation	Knowledge about health hazard is incomplete.	

Innalation	Knowledge about neath nazard is incomplete.
Skin contact	Knowledge about health hazard is incomplete.
Eye contact	Causes eye irritation. Endocrine effects.
Ingestion	Based on information from therapeutic use, this material may cause:
Symptoms related to the physical, chemical, and toxicological characteristics	Corticosteroids: Fluid and electrolyte imbalance. Adrenal suppression. Immunosuppression. Cushing's syndrome. High blood pressure. Gastrointestinal disturbances. Headache. Lightheadedness. Weakness. Visual disturbances. Mood or mental status changes. Infection. Thinning skin. Swelling. Bruising. Bone fractures. Back pain. Joint pain. Tremors. Menstrual irregularities. Impotence.

Information on toxicological effects

Acute toxicity

Acute toxicity	Species		Toot Dooulto
Components	Species		Test Results
Dibasic Calcium Phosphate (CAS	5 7757-93-9)		
<u>Acute</u>			
Dermal LD50	Dabbit		> 7040 mg/kg
LD30	Rabbit		> 7940 mg/kg
			> 2000 mg/kg, 24 Hours
Inhalation			
Dust	Det		
LC50	Rat		> 2.6 mg/l, 4 Hours
Oral			
LD50	Rat		> 10000 mg/kg
			> 2000 mg/kg
Magnesium Stearate (CAS 557-0	4-0)		
<u>Acute</u>			
Inhalation			
LC50	Rat		> 2 mg/l
Microcrystalline Cellulose (CAS 9	004-34-6)		
Acute			
Dermal			
LD50	Rabbit		> 2 g/kg
Inhalation			
LC50	Rat		> 5.05 mg/l, 4 hours
Oral			
LD50	Rat		> 5 g/kg
Stearic Acid (CAS 57-11-4)			
Oral			
LD50	Rat		4.6 g/kg
<u>Acute</u>			
Dermal			
LD50	Rabbit		> 5000 mg/kg
Skin corrosion/irritation	Knowledge about health haz	ard is incomplete.	
Serious eye damage/eye irritation	Causes eye irritation.	·	
Local effects			
Magnesium Stearate		Eye irritation Result: Negative.	
Microcrystalline Cellulose		Species: Rabbit Eye irritation	
Changin Anid		Result: Negative. Species: Rabbit	
Stearic Acid		Eye irritation Result: Negative. Species: Rabbit	
Dibasic Calcium Phosphate		Irritancy test Result: Mild Species: Rabbit Organ: Eye Irritancy test Result: Non-irritant Species: Rabbit Organ: Skin	

	Local effects		
	Stearic Acid		Skin irritation
			Result: Mild. Species: Human
			Test Duration: 3 day
	Magnesium Stearate		Skin irritation
			Result: Negative.
	Microcrystalline Cellulose		Species: Rabbit Skin irritation
			Result: Negative.
			Species: Rabbit
	Stearic Acid		Skin irritation Result: Negative.
			Species: Rabbit
			Test Duration: 24 hours
_			Observation Period: 72 hours
Res	spiratory or skin sensitization		
	Respiratory sensitization	Knowledge about health hazar	
	Skin sensitization	Knowledge about health hazar	d is incomplete.
	Stearic Acid		7 % Skin sensitization
			Result: Negative. Species: Human
	Microcrystalline Cellulose		Skin sensitization
	-		Result: Negative.
-			Species: Guinea pig
Gei	m cell mutagenicity	Knowledge about mutagenicity	is incomplete.
	Mutagenicity Microcrystalline Cellulose		Ames test
	Microcrystalline Cellulose		Result: Negative.
	Stearic Acid		Ames test
	Dradajagaa		Result: Negative.
	Prednisone		Ames test (Salmonella typhimurium) Result: Negative
	Microcrystalline Cellulose		Micronucleus test
	Dradniaana		Result: Negative.
	Prednisone		Mouse lymphoma assay Result: Negative
	Microcrystalline Cellulose		Mouse lymphoma assay
	Dradajagaa		Result: Negative.
	Prednisone		Mutagenicity, In vivo chromosome aberration studies in rat bone marrow
			Result: Negative
	Stearic Acid		Mutagenicity: Induction of crossovers and aneuploidy in S.
			cerevisiae. Result: Negative
	Microcrystalline Cellulose		Unscheduled DNA synthesis
			Result: Negative.
Car	cinogenicity	Knowledge about carcinogenic	· ·
	Stearic Acid		0.3 % Carcinogenicity Result: Negative.
			Species: Rat
			Test Duration: 209 days
	Microcrystalline Cellulose		30 % Oral feed study Result: Negative.
			Species: Rat
			Test Duration: 72 weeks
	Prednisone		5 mg/kg/day Carcinogenicity Result: Negative
			Species: Mouse
	Stearic Acid		50 g/kg Oral feed study
			Result: Negative. Species: Mouse
	Magnesium Stearate		Carcinogenicity: Implanted mouse bladders
	v		Result: Negative.

Carcinogenicity		
Microcrystalline Cellulose	Carcinogeni Result: Neg Species: Ra Test Duratio	t
IARC Monographs. Overall E	valuation of Carcinogenicity	
Prednisone (CAS 53-03-2) 3 Not classi	fiable as to carcinogenicity to humans.
OSHA Specifically Regulated	l Substances (29 CFR 1910.1001-1050)	
Not regulated.		
US. National Toxicology Pro	gram (NTP) Report on Carcinogens	
Not listed.		
Reproductive toxicity	Suspected of damaging fertility or the unborn child. Most studies have concluded that therapeutic use of corticosteroids by pregnant women does not cause adverse effects on the fetus. A small increase in the incidence of cleft palate was seen in some human studies. Infants born to mothers who received substantial doses of corticosteroids during pregnancy should be observed for signs of hypoadrenalism.	
Reproductivity		
Prednisone		
Microcrystalline Cellulose		ay Reproductivity / developmental adverse effects on the offspring observed. t
Specific target organ toxicity - single exposure	Knowledge about health hazard is incomplete.	
Specific target organ toxicity - repeated exposure	Causes damage to organs (endocrine system) through prolonged or repeated exposure.	
Aspiration hazard	Based on available data, the classification criteria are not met.	
Further information	Potent pharmacologically active material. Occupational exposure to small amounts may cause physiological effects.	

12. Ecological information

Components		Species	Test Results
Dibasic Calcium Phosphate	(CAS 7757-9	93-9)	
Aquatic			
Acute			
Crustacea	EC50	Daphnia magna	> 1000 mg/l
Prednisone (CAS 53-03-2)			
Aquatic			
Acute			
Algae	IC50	Algae	31 mg/l, 72 hours
Stearic Acid (CAS 57-11-4)			
Acute			
Other	EC50	Pseudomonas putida	> 100 mg/l, 16 hours
Aquatic			
Acute			
Fish	LC50	Carp (Cyprinus carpio)	> 1000 mg/l, 48 hours
rsistence and degradability	No data is available on the degradability of this product.		
paccumulative potential			
Octanol/water partition co	efficient log	Kow	
Prednisone	-	1.46	
bility in soil	No data available.		
her adverse effects	No other adverse environmental effects (e.g. ozone depletion, photochemical ozone creation potential, endocrine disruption, global warming potential) are expected from this component.		

13. Disposal considerations

Disposal instructions	Dispose in accordance with all applicable regulations. Under RCRA, it is the responsibility of the user of the product to determine, at the time of disposal, whether the product meets RCRA criteria for hazardous waste.
Local disposal regulations	Dispose in accordance with all applicable regulations.
Hazardous waste code	The waste code should be assigned in discussion between the user, the producer and the waste disposal company.
Waste from residues / unused products	Empty containers or liners may retain some product residues. This material and its container must be disposed of in a safe manner (see: Disposal instructions).
Contaminated packaging	Empty containers should be taken to an approved waste handling site for recycling or disposal. Since emptied containers may retain product residue, follow label warnings even after container is emptied.

14. Transport information

DOT

Not regulated as dangerous goods.

ΙΑΤΑ

Not regulated as dangerous goods.

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	Not applicable.
General information	It is the shipper's responsibility to determine the correct transport classification at the time of shipment.

15. Regulatory information

US federal regulations

One or more components are not listed on TSCA.

CERCLA/SARA Hazardous Substances - Not applicable. This product is a "Hazardous Chemical" as defined by the OSHA Hazard Communication Standard, 29 CFR 1910.1200.

TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D)

Not regulated.

CERCLA Hazardous Substance List (40 CFR 302.4)

Not listed.

SARA 304 Emergency release notification

Not regulated.

OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

Not regulated.

Superfund Amendments and Reauthorization Act of 1986 (SARA)

Hazard categories	Immediate Hazard - Yes Delayed Hazard - Yes
	Fire Hazard - No
	Pressure Hazard - No
	Reactivity Hazard - No

SARA 302 Extremely hazardous substance

Not listed.

SARA 311/312 Hazardous Yes chemical

SARA 313 (TRI reporting) Not regulated.

Other federal regulations

Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List

Not regulated.

Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130)

Not regulated.

Safe Drinking Water Act Not regulated. (SDWA)

US state regulations

California Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65): This material is not known to contain any chemicals currently listed as carcinogens or reproductive toxins.

International Inventories

Country(s) or region	Inventory name	On inventory (yes/no)*
Australia	Australian Inventory of Chemical Substances (AICS)	Yes
Canada	Domestic Substances List (DSL)	No
Canada	Non-Domestic Substances List (NDSL)	Yes
China	Inventory of Existing Chemical Substances in China (IECSC)	No
Europe	European Inventory of Existing Commercial Chemical Substances (EINECS)	No
Europe	European List of Notified Chemical Substances (ELINCS)	No
Japan	Inventory of Existing and New Chemical Substances (ENCS)	Yes
Korea	Existing Chemicals List (ECL)	No
New Zealand	New Zealand Inventory	No
Philippines	Philippine Inventory of Chemicals and Chemical Substances (PICCS)	Yes
United States & Puerto Rico	Toxic Substances Control Act (TSCA) Inventory	Yes

*A "Yes" indicates that all components of this product comply with the inventory requirements administered by the governing country(s) A "No" indicates that one or more components of the product are not listed or exempt from listing on the inventory administered by the governing country(s).

16. Other information, including date of preparation or last revision

Issue date Revision date	09-08-2010 11-13-2018
Version #	04
Disclaimer	USP Reference Standards are sold for chemical test and assay purposes only, and NOT for human consumption. The information contained herein is applicable solely to the chemical substance when used as a USP Reference Standard and does not necessarily relate to any other use of the substance described, (i.e. at different concentrations, in drug dosage forms, or in bulk quantities). USP Reference Standards are intended for use by persons having technical skill and at their own discretion and risk. This information has been developed by USP staff from sources considered reliable but has not been independently verified by the USP. Therefore, the USP Convention cannot guarantee the accuracy of the information in these sources nor should the statements contained herein be considered an official expression. NO REPRESENTATION OR WARRANTY, EXPRESS OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE is made with respect to the information contained herein.